Monoclonal Mouse Anti-Human CD246, ALK Protein
Clone ALK1
Code No./ Code/ Code-Nr. M 7195
Edition/ Ausgabe 19.12.02

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD246, ALK Protein, Clone ALK1, is intended for use in immunocytochemistry. The antibody labels normal human ALK protein and the NPM-ALK chimeric protein, and is a useful tool for the identification of the subgroup of anaplastic large-cell lymphomas (ALCL) that are ALK positive (1, 2). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	CD246 (3).
Introduction	The wild-type anaplastic lymphoma kinase (ALK) protein is a 200 kDa transmembrane receptor tyrosine kinase (4). Postnatal ALK expression is restricted to a few scattered cells in the nervous system (some glial cells and neurons, and a few endothelial cells and pericytes) (1). About 72.5% of ALK-positive ALCL are associated with a (2;5) chromosomal translocation, where the nucleophosmin (NPM) gene located at 5q35 fuses with the ALK gene located at 2p23. With the fusion, the portion of the NPM gene encoding the N-terminal part of the NPM protein is juxtaposed to the part of the ALK gene that codes for the entire cytoplasmic region of the ALK protein. As a consequence, the ALK gene comes under the control of the NPM promoter, which induces a permanent and ubiquitous transcription of the NPM-ALK hybrid gene, resulting in the production of a 80 kDa NPM-ALK chimeric protein (4). In 3 large series of ALCL, 15-28% of chimeric ALK-positive lymphomas were negative for the t(2;5) translocation, and the main alternative fusion gene was identified as the tropomyosin 3 gene representing 17.5% of the ALK-positive ALCL cases (4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> ALK1 (1). <u>Isotype:</u> IgG3, kappa. <u>Mouse IgG concentration:</u> see label on vial.
Immunogen	DHFR-ALK recombinant protein comprising mouse dihydrofolate reductase (DHFR) and amino acids 1359-1460 of the full length human ALK protein, corresponding to amino acids 419-520 of the chimeric NPM-ALK protein (1).
Specificity	The antibody was clustered as anti-CD246 at the Seventh International Workshop and Conference on Human Leucocyte Differentiation Antigens held in Harrogate in 2000 (3). In Western blotting of lysate of the t(2;5)+ SU-DHL-1 cell line, the antibody labels a band of 80 kDa corresponding to the NPM-ALK protein. The antibody also recognizes a band of 200 kDa corresponding to the ALK protein in the human rhabdomyosarcoma cell line Rh30, that expresses full length ALK protein, and in the 293T cell line transfected with cDNA encoding full length ALK protein. No reactivity was observed when the 293T cell line was transfected with the vector alone (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or B5-fixative (4). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, DakoCytomation Target Retrieval Solution, code No. S 1700, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections and cell preparations (5).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD246, ALK Protein, code No. M 7195, may be used at a dilution range of 1:25-1:50 when applied on formalin-fixed, paraffin-embedded sections of t(2;5)+ ALCL and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval solution, High pH, code No. S 3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic and/or nuclear staining pattern.
(102304-001)	M 7195/EFG/CE/19.12.02 p. 1/4

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
- Pulford K, Lamant L, Morris SW, Butler LH, Wood KM, Stroud D, et al. Detection of anaplastic lymphoma kinase (ALK) and nuclear protein nucleophosmin (NPM)-ALK proteins in normal and neoplastic cells with the monoclonal antibody ALK1. Blood 1997;89:1394-404.
- Gascoyne RD, Aoun P, WU D, Chhanabhai M, Skinnider BF, Greiner TC, et al. Prognostic significance of anaplastic lymphoma kinase (ALK) protein expression in adults with anaplastic large cell lymphoma. Blood 1999;93:3913-21.
- Pulford K, Mason DY. TC5. CD246 (ALK) cluster report. In: Mason D, André P, Bensussan A, Buckley C, Civin C, Clark E, et al., editors. Leucocyte typing VII. White cell differentiation antigens. Proceedings of the 7th International Workshop and Conference; 2000 Jun 19-23; Harrogate, United Kingdom. New York: Oxford University Press Inc.; 2002. p. 694-6.
- Stein H, Foss H-D, Dürkop H, Marafioti T, Delsol G, Pulford K, et al. CD30⁺ anaplastic large cell lymphoma: a review of its histopathologic, genetic, and clinical features. Blood 2000;96:3681-95.
- Lamant L, Pulford K, Bischof D, Morris SW, Mason DY, Delsol G, et al. Expression of the ALK tyrosine kinase gene in neuroblastoma. Am J Pathol 2000;156:1711-21.
- Delsol G, Lamant L, Mariamé B, Pulford K, Dastugue N, Brousset P, et al. A new subtype of large B-cell lymphoma expressing the ALK kinase and lacking the 2;5 translocation. Blood 1997;89:1483-90.

Explanation of symbols/ Légende des symboles/ Erläuterung der Symbole			
REF	Catalogue number Référence du catalogue Bestellnummer	 2 °C – 8 °C Temperature limitation Limites de température Zulässiger Temperaturbereich	 Manufacturer Fabricant Hersteller
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-Vitro-Diagnostikum	LOT Batch code Code du Lot Chargenbezeichnung	
	Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Use by Utiliser jusque Verwendbar bis	

(102304-001)	M 7195/EFG/CE/19.12.02 p. 4/4	
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17		