**Monoclonal Mouse  
Anti-Human  
Androgen Receptor  
Clone AR441**

**ENGLISH**  
**Code M3562**

**Intended use**

For In Vitro Diagnostic Use.

**Summary and explanation**

Androgen receptor is an intracellular protein which belongs to a large family of hormone-induced transcription factors.<sup>2,3</sup> Androgen receptor mediates the biological actions of physiological androgens such as testosterone and 5  $\alpha$ -dihydrotestosterone, which are essential for differentiation, development and maintenance of the male reproductive organs. Complexes of androgen and receptor regulate growth responses in accessory sex organs by modulating specific gene transcription.<sup>3</sup>

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Anti-androgen receptor is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2, and 0.015 mol/L sodium azide.

Clone: AR441      Isotype: IgG<sub>1</sub>, kappa  
Mouse IgG concentration mg/L: See label on vial.

Anti-human Androgen Receptor may be used at a dilution of 1:50 to 1:75 in the LSAB™+ method determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

**Immunogen**

Synthetic peptide corresponding to amino acids 229–315 of the human androgen receptor coupled to keyhole limpet hemocyanin<sup>1</sup>

**Specificity**

Monoclonal mouse anti-human androgen receptor, clone AR441 (anti-androgen receptor) specifically immunostains the nuclei of positive cells in human tissues. In Western blot analysis, the antibody identifies a 110 and 112 kD doublet in extracts of the metastatic prostate cancer cell line LNCap and in extracts of cells transfected with the gene for androgen receptor.<sup>1</sup> The antibody is unreactive with mouse androgen receptor in immunoblots of prostate extracts. Anti-androgen receptor is capable of recognizing and forming complexes with native androgen bound receptor.<sup>1</sup>

**Materials required, but not supplied**

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

**Storage**

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

**Specimen preparation**

*Paraffin Sections*

Anti-androgen receptor can be used on formalin-fixed, paraffin-embedded tissue sections.

The deparaffinized tissue sections must be treated with heat prior to the IHC staining procedure. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended.

When using the water bath method, preheat a Coplin jar containing Target Retrieval Solution, High pH (code S3307 or S3308), as well as a water bath to 95–99 °C. When the temperature has stabilized, place tissue sections in the Coplin jar containing the preheated buffer. Heat the tissue sections for 20 minutes. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with distilled water and place slides in buffer.

### *Cryostat Sections and Cell Smears*

Anti-androgen receptor can be used for labeling acetone-fixed cryostat sections or fixed cell smears.

### **Staining procedure**

Follow the recommended procedure for the detection system selected.

### **Staining interpretation**

The cellular staining pattern for anti-androgen is nuclear.

### **Performance characteristics**

#### *Normal Cells*

Androgen receptor has been localized by immunohistochemistry (IHC) to the nuclei of a variety of normal cells. In the prostate, nuclear immunoreactivity has been observed in the secretory cells, basal cells, smooth muscle cells and fibroblasts. Moderate immunostaining has been found in spermatogonia, spermatocytes and Sertoli cells of the testes, with weaker staining seen in Leydig and peritubular myoid cells. Strong to moderate expression of androgen receptor has been demonstrated by IHC in uterine myometrial cells, cervical squamous cells, ovarian stromal cells, and syncytiotrophoblasts and cytotrophoblasts of the placenta.<sup>4</sup> Endometrial glandular and stromal cells also express androgen receptor at all phases of the menstrual cycle.<sup>4,5</sup> In the breast, the nuclei of acinar cells, myoepithelial cells and stromal cells were reported to stain positively. Androgen receptor was also demonstrated in cardiac muscle, striated muscle and in smooth muscle of arteries, bladder and GI tract. In skin, positive staining has been observed in cells of the hair follicles, sweat glands and squamous cells of the epidermis. Other cells which have been shown to give moderate to strong immunostaining include hepatocytes of the liver and cells of the anterior and posterior pituitary.<sup>4</sup>

#### *Tumor Cells*

In adenocarcinomas of the prostate androgen receptor has been demonstrated by IHC to be present in the nuclei of most neoplastic cells with neuroendocrine differentiation.<sup>6</sup> The majority of breast and ovarian carcinomas evaluated also express androgen receptor, as determined by IHC and biochemical assays.<sup>7-10</sup> The nuclei of glandular and stromal endometrial cells have been shown to stain positively in adenomyosis and external endometriosis. Malignant cells in both glandular and solid lesions of grade II endometrial adenocarcinomas were found to be immunoreactive.<sup>5</sup> Androgen receptor expression has also been demonstrated in meningiomas by Northern blot analysis and immunohistochemistry.<sup>11</sup>

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### Références

#### Literatur

1. Report on file.
2. Zhuang YH, et al. Subcellular location of androgen receptor in rat prostate seminal vesicle and human osteosarcoma MG-63 cells J Steroid Biochem Mol Biol 1992; 41(3-8):693
3. Jänne OA, et al. Androgen receptor and mechanism of androgen action. Ann Med 1993; 25:83
4. Kimura N, et al. Immunocytochemical localization of androgen receptor with polyclonal antibody in paraffin-embedded human tissues. J Histochem Cytochem 1993; 41(5):671
5. Horie K, et al. Immunohistochemical localization of androgen receptor in the human endometrium, decidua, placenta and pathological conditions of the endometrium. Hum Reproduct 1992; 7(10):1461
6. Nakada SY, et al. The androgen receptor status of neuroendocrine cells in human benign and malignant prostatic tissue. Canc Res 1993; 53:1967
7. Schulz TJ, et al. Immunocytochemical analysis of androgen receptor in breast cancer and correlation with estrogen and progesterone receptors. Lab Invest 1993; 68(1):19A(#92)
8. Bryan RM, et al. Androgen receptors in breast cancer. Canc 1984; 54(11):2436
9. Chadha S, et al. An immunohistochemical evaluation of androgen and progesterone receptors in ovarian tumors. Hum Pathol 1993; 24(1):90
10. Kühnel R, et al. Androgen receptor predominance in human ovarian carcinoma. J Steroid Biochem 1987; 26(3):393
11. Carroll RS, et al. Androgen receptor expression in meningiomas. J Neurosurg 1995; 82:453

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

  
Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA

 Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark

Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

Tel +45 4485 9500  
Fax +45 4485 9595

www.dako.com

PT0039/Rev C