ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human BCL2 Oncoprotein, Clone 124, is intended for use in immunohistochemistry. The antibody labels cells expressing BCL2 oncoprotein. Positive results aid in the classification of follicular lymphomas and various diffuse lymphoproliferative diseases. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	BCL2 oncoprotein is a blocker of apoptotic cell death. Gene transfer experiments have shown that elevated levels of this protein can protect a wide variety of cells from diverse cell death stimuli ranging from growth factor withdrawal and cytotoxic lymphokines to virus infection and DNA-damaging, anticancer drugs and radiation (2, 3). BCL2 oncoprotein resides on the cytoplasmic side of the mitochondrial outer membrane, endoplasmic reticulum and nuclear envelope (2, 4), and has a molecular mass of 26 kDa (3). The <i>BCL2</i> gene is involved in the t(14;18) chromosomal translocation found in 85% of human follicular lymphomas and 20% of diffuse B-cell lymphomas (4). In this translocation, the <i>BCL2</i> gene at chromosome segment 18q21 is juxtaposed with the Ig heavy chain locus at 14q32, resulting in deregulated expression of BCL2 oncoprotein (4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 124 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is adjusted to ascertain comparable lot-to-lot immunohistochemical staining performance with a reference lot.
Immunogen	Synthetic peptide comprising amino acids 41-54 of human BCL2 oncoprotein (1).
Specificity	In Western blotting of extracts of normal human spleen (1, 5), t(14;18)-positive follicular lymphoma (1), and myloid leukaemic cell lines (5) the antibody labels solely a band of 26 kDa, corresponding to BCL2 oncoprotein under both nonreducing (1) and reducing conditions (1, 5). The antibody labels the myeloid leukaemic cell lines, HL-60 (promyelocytic), KG1 (myeloblastic), GM-1 (monoblastic) and K562 (erythromyeloid) (5).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin or Bouin's (5). Pre-treatment of tissues with heat-induced epitope retrieval is required. For tissues fixed in formalin, optimal results are obtained with 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Less optimal results are obtained with Dako Target Retrieval Solution, High pH, code S3308, Dako Target Retrieval Solution, code S1700, or 10 mmol/L citrate buffer, pH 6.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections and cell smears (1).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human BCL2 Oncoprotein, code M0887, may be used at a dilution range of 1:50-1:100 when applied on formalin-fixed, paraffin-embedded sections of tonsil tissue and using 20 minutes heat-induced epitope retrieval in 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> Dako LSAB™+/HRP kit, code K0679, and Dako EnVision™+/HRP kits, codes K4004 and K4006 are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer.
Performance characteristics	Cells labelled by the antibody display a cytoplasmic staining pattern. <u>Normal tissues:</u> The antibody labels almost all peripheral blood lymphocytes. In lymphoid tissue, small lymphocytes in the mantle zones and T-cell areas are positive whereas very few cells in germinal centres are labelled. In the spleen, many cells in both T- and B-cell areas and the red pulp are labelled by the antibody. In the thymus, many cells in the medulla are labelled, while most cells in the cortex show weak or negative staining (1).

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
- Pezzella F, Tse AGD, Cordell JL, Pulford KAF, Gatter KC, Mason DY. Expression of the *bcl-2* oncogene protein is not specific for the 14;18 chromosomal translocation. Am J Pathol 1990;137:225-32.
- Adams JM, Cory S. The bcl-2 protein family: Arbiters of cell survival. Science 1998;281:1322-6.
- Kusenda J. Bcl-2 family proteins and leukaemia [minireview]. Neoplasma 1998;45:117-22.
- Yang E, Korsmeyer SJ. Molecular thanatopsis: A discourse on the BCL2 family and cell death [review]. Blood 1996;88:386-401.
- Delia D, Aiello A, Soligo D, Fontanella E, Melani C, Pezzella F, et al. bcl-2 proto-oncogene expression in normal and neoplastic human myeloid cells. Blood 1992;79:1291-8.
- Hirakawa N, Naka T, Yamamoto I, Fukuda T, Tsuneyoshi M. Overexpression of bcl-2 protein in synovial sarcoma: A comparative study of other soft tissue spindle cell sarcomas and an additional analysis by fluorescence in situ hybridization. Hum Pathol 1996;27:1060-5.
- Soini Y, Pääkkö P. bcl-2 is preferentially expressed in tumours of muscle origin but is not related to p53 expression. Histopathology 1996;28:141-5.

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