## Monoclonal Mouse Anti-Human BCL6 Protein

Clone PG-B6p  
Code M7211

ENGLISH
<p><b>Intended use</b> For in vitro diagnostic use.</p>

Monoclonal Mouse Anti-Human BCL6 Protein, Clone PG-B6p, is intended for use in immunohistochemistry. The antibody labels cells expressing BCL6 protein and is a useful tool for differentiating proliferating centres in B-cell chronic lymphocytic lymphomas (BCL2+/BCL6-) from trapped germinal centres in mantle cell lymphomas (BCL2-/BCL6+), and for identifying neoplastic cells in cases of nodular, lymphocyte-predominance Hodgkin's disease (1). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and explanation

The proto-oncogene *BCL6* encodes a 92-98 kDa POZ/zinc finger protein. Through its six C-terminal zinc-finger motifs, the BCL6 protein binds to specific DNA sequences of target genes, where it acts as a transcriptional repressor (2, 3).

*BCL6* expression may be ubiquitous, as low levels are detected in many tissues. However, high levels of *BCL6* expression are limited to cells of the B-lymphoid lineage and skeletal muscle. In the B-lymphoid cell lineage, *BCL6* expression seems to be tightly regulated during differentiation, and strong expression is observed preferentially in germinal centre (GC) B cells, but not in plasma cells, and precursor, mantle zone, and memory B cells. Neoplastic counterparts of GC B cells, including Burkitt's, follicular, and diffuse large cell lymphomas, also display high levels of *BCL6* expression (3).

*BCL6* may be implicated in chromosomal translocations where the regulatory region of the *BCL6* gene is replaced by a heterologous reciprocal partner such as the immunoglobulin (IG) genes. Promotor substitution leads to deregulation of the *BCL6* expression, which may be associated with lymphomagenesis (3). *BCL6* rearrangement is one of the most common genetic abnormalities in B-cell non-Hodgkin's lymphoma, with an especially high frequency (28-40%) in diffuse large cell lymphoma (2).

### Reagent provided

Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN<sub>3</sub>.

**Clone:** PG-B6p (1). **Isotype:** IgG1, kappa.

**Mouse IgG concentration:** see label on vial.

The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

### Immunogen

Recombinant glutathione S-transferase (GST)-BCL6 protein corresponding to amino acids 3-484 of the human BCL6 protein (1).

### Specificity

SDS-PAGE analysis of immunoprecipitates formed between the antibody and total lysates of Bjab cells (BCL6+) and BCL6 transfected EB3 cells, shows reaction with a ~95 kDa fragment corresponding to BCL6. No reaction is observed with Rd cells (BCL6-) and EB3 cells transfected with a control plasmid (1).

In Western blotting of total lysates of Bjab cells and BCL6 transfected EB3 cells, or control transfected EB3 cells, the antibody did not label any bands, suggesting lack of reactivity with denatured BCL6 protein (1).

As demonstrated by immunocytochemistry, the antibody cross-reacts with the BCL6-equivalent protein in various mammals, including cow, rabbit, rat, sheep, and swine. No cross-reactivity was observed in chicken and pigeon (1).

### Precautions

- For professional users.
- This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
- As with any product derived from biological sources, proper handling procedures should be used.
- Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
- Unused solution should be disposed of according to local, State and Federal regulations.

### Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.

### Specimen preparation

**Paraffin sections:** The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, High pH, code S3308, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. However, Dako Target Retrieval Solution, code S1700, or 10 mmol/L citrate buffer, pH 6.0 and pre-treatment of tissues with proteinase K were found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure.

**Frozen sections and cell preparations:** The antibody can be used for labelling acetone fixed, frozen sections.

### Staining procedure

**Dilution:** Monoclonal Mouse Anti-Human BCL6 Protein, code M7211, may be used at a dilution range of 1:10-1:20 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, High pH, code S3308, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen.

**Visualization:** Dako LSAB™+/HRP kit, code K0679, and Dako EnVision™+/HRP kits, code K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

**Automation:** The antibody is well-suited for immunohistochemical staining using automated platforms.

### Performance characteristics

Cells labelled by the antibody show a diffuse/microgranular staining confined to the nucleus. Intense cytoplasmic staining unassociated with metaphase chromosomes is noted only in mitotic figures (1).

**Normal tissues:** In tonsil, the antibody strongly labels the nuclei of centroblasts in the dark zone and centrocytes in the basal and apical light zones of germinal centres. In addition the antibody labels approximately 10% of germinal centre T cells. No labelling was observed in plasma cells, macrophages, follicular dendritic cells, and IgD+/IgM+ follicular mantle lymphocytes. In spleen, the staining pattern was similar to that of tonsil, but in addition the antibody labelled a few scattered lymphoid-like cells of undefined phenotype. In extra-lymphoid tissues, the antibody faintly labelled squamous epithelia in the tonsil, thymus, and skin. No labelling was observed in liver, thyroid, striated muscle, or kidney (1).

**Abnormal tissues:** 173 cases of human lymphoid neoplasms were tested with the antibody. Of B-cell lymphomas or leukaemias, 24/24 follicle centre lymphomas, 29/30 diffuse large cell lymphomas, and 13/13 Burkitt's lymphomas, were positive. No labelling was observed in 7/7 pre-B acute lymphoblastic leukaemias, 16/16 small lymphocytic/B-CLL, 6/6 hairy cell leukaemias, 22/22 mantle cell leukaemias, and 14/14 marginal zone lymphomas. Of T-cell lymphomas or leukaemias, 4/8 anaplastic large cell lymphomas were positive, whereas no labelling was observed in 10/10 acute lymphoblastic lymphomas, 3/3 mycosis fungoides, 6/6 peripheral T-cell lymphomas, and 2/2 peripheral T-cell (AILD-like) lymphomas. In Hodgkin's disease, 5/5 nodular, lymphocyte predominance, 1/4 nodular sclerosis, and 1/3 mixed cellularity types were positive (1).

In systemic AIDS-related non-Hodgkin's lymphoma, the antibody labelled 11/11 small noncleaved cell lymphomas, 7/7 large noncleaved cell lymphomas (LNCCL), and 2/9 large cell immunoblastic lymphoma plasmacytoid (IBLP). In AIDS-primary central nervous system lymphomas, the antibody labelled 4/4 LNCCL and none of 4 cases of IBLP. No labelling was observed in 5 cases of primary effusion lymphoma (PEL) and in 6 AIDS-PEL cell lines (4)

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