**Monoclonal Mouse Anti-Human CD138**  
Clone MI15  
**Code No./ Code/ Code-Nr. M 7228**  
Edition/ Ausgabe 01.03.04

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human CD138, Clone MI15, is intended for use in immunocytochemistry together with a panel of other antibodies. CD138 is present in malignancies involving terminally differentiated plasma cells, and it is associated with multiple myeloma (1, 2, 3, 4). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.
<b>Synonym for antigen</b>	Syndecan-1 (1).
<b>Introduction</b>	CD138 is a transmembrane proteoglycan, bearing heparan sulfates and sometimes chondroitin sulfates. It consists of a 30.5 kDa core protein and five glycosaminoglycans. CD138 is a member of the syndecan family. The reported Mr of CD138 ranges from 85 000 to 92 000, presumably depending on the glycosaminoglycan pattern of the molecule (1).  Within the haematopoietic system, CD138 is mainly confined to late stages of B-cell differentiation. It is typically expressed at high cellular density in normal and malignant plasma cells (5, 6) and in lymphoplasmacytoid cells (6). In mature tissue CD138 is expressed by simple and stratified epithelia, fibroblasts, stratified keratinocytes and endothelial cells (1).  The function of CD138 is largely unknown. However, the intracellular tail of CD138 seems to bind together the cytoskeleton with cellular components (1), and the extracellular part of the molecule seems to bind different ligands. Presumably, CD138 plays a role in cellular functions, such as proliferation, programmed cell death, cell-matrix and cell-cell adhesion (5), for example it mediates adhesion of myeloma cells to type I collagen (2).  CD138 expression is reduced during malignant transformation of various epithelia (7), and the antigen is rapidly lost (shed) by myeloma cells entering into apoptosis, thus CD138 is a marker of viable myeloma cells (1, 2). CD138 is lost upon malignization of several stratified epithelia (7). Together with BCL-6 and IRF4, CD138 may be of value for the diagnosis of AIDS-related non-Hodgkin lymphomas and HIV-related Hodgkin lymphomas (8).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> .  <u>Clone:</u> MI15. <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
<b>Immunogen</b>	A mixture of U266 and XG-1 human myeloma cell lines (3, 9).
<b>Specificity</b>	Anti-CD138, MI15, was included in the Sixth International Workshop and Conference on Human Leucocyte Differentiation Antigens (Kobe 1996), and studies by a number of different laboratories confirmed its reactivity with CD138 (1).  The epitope recognized by anti-CD138, MI15, is found within the ectodomain (extracellular part) of the CD138 core protein. Blocking experiments have revealed that the epitope recognized by MI15 seems to be partly overlapping the epitope recognized by another antibody to CD138, clone B-B4 (1, 5). Anti-CD138, MI15, reacts strongly with CD138 expressed by multiple myeloma-derived cell lines, such as U266 and XG-1 (9).
<b>Precautions</b>	1. For professional users.  2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.  3. As with any product derived from biological sources, proper handling procedures should be used.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used on formalin-fixed, paraffin-embedded tissue sections. Antigen retrieval by heating in DakoCytomation Target Retrieval Solution, code No. S 1700 or S 1699 is mandatory. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen tissue sections, fixed cell smears, and fixed cytopins.  The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD138, code No. M 7228, may be used at a dilution range of 1:25-1:50 when applied on paraffin sections, acetone-fixed frozen tissue sections, fixed cell smears, and fixed cytopins, using 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931 diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No. K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.
<b>Limitations</b>	Owing to the sensitivity of CD138 towards proteolytic cleavage, the processing of samples, such as freezing or storage before staining with the antibody, may play an important role in the outcome of the staining (1).
(102471-002)	M 7228/EFG/KIM/01.03.04 p. 1/4
DakoCytomation Denmark A/S · Produktionsvej 42 · DK-2600 Glostrup · Denmark · Tel. +45 44 85 95 00 · Fax +45 44 85 95 95 · CVR No. 33 21 13 17	

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
- Wijdenes J, Clément C, Klein B, Dore J-M. New B-cell CD antigens. BC29: CD138 (syndecan-1) workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEG, Goyert SM, Mason DY, Miyasaka M, et al., editors. Leucocyte typing VI. White cell differentiation antigens. Proceedings of the 6th International Workshop and Conference; 1996 Nov 10-14; Kobe, Japan. New York, London: Garland Publishing Inc.; 1997. p. 249-52.
- Jourdan M, Ferlin M, Legouffe E, Horvathova M, Liautard J, Rossi JF, et al. The myeloma cell antigen syndecan-1 is lost by apoptotic myeloma cells. Br J Haematol 1998;100:637-46.
- Costes V, Magen V, Legouffe E, Durand L, Baldet P, Rossi J-F, et al. The MI15 monoclonal antibody (anti-syndecan-1) is a reliable marker for quantifying plasma cells in paraffin-embedded bone marrow biopsy specimens. Hum Pathol 1999;30:1405-11.
- van Zaanen HCT, Vet RJWM, de Jong CM, von dem Borne AEGKr, van Oers MHJ. A simple and sensitive method for determining plasma cell isotype and monoclonality in bone marrow using flow cytometry. Br J Haematol 1995;91:55-9.
- Gattei V, Godeas C, Degan M, Rossi FM, Aldinucci D, Pinto A. Characterization of anti-CD138 monoclonal antibodies as tools for investigating the molecular polymorphism of syndecan-1 in human lymphoma cells. Br J Haematol 1999;104:152-62.
- Sebestyén A, Berdzi L, Mihalik R, Paku S, Matolcsy A, Kopper L. Syndecan-1 (CD138) expression in human non-Hodgkin lymphomas. Br J Haematol 1999;104:412-9.
- Inki P, Jalkanen M. The role of syndecan-1 in malignancies [review]. Ann Med 1996;28:63-7.
- Carbone A, Gloghini A, Larocca LM, Capello F, Pierconti F, Canzonieri V, et al. Expression profile of MUM1/IRF4, BCL-6, and CD138/syndecan-1 defines novel histogenetic subsets of human immunodeficiency virus-related lymphomas. Blood 2001; 97:744-51.
- Horvathova M, Gaillard JP, Liautard J, Duperray C, Lavabre-Bertrand T, Bourquard P, et al. B30.6. Identification of novel and specific antigens of human plasma cells by mAb. In: Schlossman SF, Bousmell L, Gilks W, Harlan JM, Kishimoto T, Morimoto C, et al., editors. Leucocyte typing V. White cell differentiation antigens. Proceedings of the 5th International Workshop and Conference; 1993 Nov 3-7; Boston, USA. Oxford, New York, Tokyo: Oxford University Press; 1995. p. 713-4.

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