**Monoclonal Mouse  
Anti-Human  
Caldesmon  
Clone h-CD**

---

**ENGLISH**  
**Code M3557**

**Intended use**

For In Vitro Diagnostic Use.

**Summary and explanation**

Caldesmon is a developmentally regulated protein involved in smooth muscle and non-muscle contraction.<sup>2,3</sup>

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

**Reagent provided**

Anti-human caldesmon, h-CD is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide.

Clone: h-CD<sup>1</sup> Isotype: IgG<sub>1</sub>, kappa  
Mouse IgG concentration mg/L: See label on vial.

Anti-caldesmon, h-CD may be used at a dilution of 1:50 to 1:100 in the LSAB method determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

**Immunogen**

Crude human uterus extract<sup>1</sup>

**Specificity**

Two closely related variants of human caldesmon have been identified which differ in their electrophoretic mobility and cellular distribution. The *h*-caldesmon variant (120–150 kD) is predominantly expressed in smooth muscle whereas *l*-caldesmon (70–80 kD) is found in non-muscle tissue and cells. Neither of the two variants have been detected in skeletal muscle.<sup>2</sup> Monoclonal anti-caldesmon, h-CD, recognizes only the 150 kD variant (*h*-caldesmon) in Western blots of human aortic media extracts and is unreactive with fibroblast extracts from cultivated human foreskin.<sup>1</sup>

**Materials required, but not supplied**

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

**Storage**

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

**Specimen preparation**

*Paraffin Sections*

Anti-caldesmon, h-CD can be used on formalin-fixed, paraffin-embedded tissue sections.

Prior to the IHC staining procedure, the deparaffinized tissue sections must be treated with a proteolytic enzyme followed by target retrieval. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Deparaffinized tissue sections must first be treated for 5 to 10 minutes with a mild enzyme solution. A recommended proteolytic enzyme is Proteinase K (code S3004) which must be further diluted 1:500 in 0.05 mol/L Tris-HCl, pH 7.6 to give a final concentration of 0.04 mg/mL.

Following proteolytic digestion, tissue sections must be treated with heat. When using the water bath method, preheat a Coplin jar containing 0.01 mol/L citrate buffer, pH 6.0 as well as a water bath to 95–99 °C. When the temperature has stabilized, place tissue sections into the Coplin jar containing the preheated buffer. Heat the tissue sections for 40 minutes. For improved staining results and a shorter incubation time, Target Retrieval Solution (code S1700) can be used in place of the 0.01 mol/L citrate buffer. Under these conditions the incubation time in the water bath may be reduced to 20 minutes. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with distilled water and place slides into buffer.

#### *Cryostat Sections and Cell Smears*

Anti-caldesmon, h-CD can also be used to label cryostat sections or cell smears.

#### **Staining procedure**

Follow the recommended procedure for the detection system selected.

#### **Staining interpretation**

The staining pattern for anti-caldesmon is cytoplasmic.

#### **Performance characteristics**

##### *Normal Cells*

Monoclonal anti-caldesmon was found to localize *h*-caldesmon in cryostat sections of human tissues including developing visceral smooth muscle from 10–20 week-old fetal trachea, esophagus, jejunum and uterus. Aortic smooth muscle cells from 10 and 20 week-old fetuses were unreactive with anti-caldesmon.<sup>1</sup> In the adult, caldesmon expression was demonstrated in vascular and visceral smooth muscle cells but not epithelial cells, endothelial cells or connective tissue fibroblasts. Cells of the tunica media and a subpopulation of smooth muscle cells of the sub-endothelial intima of the adult aorta were shown to stain positively with anti-caldesmon.<sup>1</sup> In immunohistochemical (IHC) studies on cryostat sections of normal human breast, monoclonal anti-caldesmon labeled the smooth muscle cells of blood vessels and a subset of myoepithelial cells in the galactophorous sinuses. Myoepithelial cells of the lobules, ducts and luminal epithelial cells of normal breast did not react positively.<sup>4</sup>

##### *Tumor Cells*

In IHC studies on breast carcinoma, monoclonal anti-caldesmon has been demonstrated to label a subpopulation of myoepithelial cells but is unreactive with myofibroblasts and tumor cells.<sup>5</sup>

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### Références

#### Literatur

1. Frid MG, et al. Phenotypic changes of human smooth muscle cells during development: Late expression of heavy caldesmon and calponin. Dev Biol 1992; 153:185
2. Sobue K and Sellers JR. Caldesmon, a novel regulatory protein in smooth muscle and nonmuscle actomyosin systems. J Biol Chem 1991; 266(19):12115
3. Glukhova MA, et al. Developmental changes in expression of contractile and cytoskeletal proteins in human aortic smooth muscle. J Biol Chem 1990; 265(22):13042
4. Lazard D, et al. Expression of smooth muscle-specific proteins in myoepithelium and stromal myofibroblasts of normal and malignant human breast tissue. Proc Natl Acad Sci USA 1993; 90:999
5. Wang NP, et al. Monoclonal antibodies (MAbs) to novel myoepithelium-associated proteins can distinguish between benign and malignant lesions of the breast. US Canad Acad Pathol Abstr 1996; 26A:135

|   |   |  |
|---|---|--|
| <b>REF</b><br>Catalogue number<br>Référence du catalogue<br>Bestellnummer   |  Temperature limitation<br>Limites de température<br>Zulässiger Temperaturbereich |  Consult instructions for use<br>Consulter les instructions d'utilisation<br>Gebrauchsanweisung beachten |
|  Manufacturer<br>Fabricant<br>Hersteller   | <b>LOT</b><br>Batch code<br>Code du lot<br>Chargenbezeichnung   |  Use by<br>Utiliser jusque<br>Verwendbar bis  |
| <b>EC REP</b><br>Authorized representative in the European Community<br>Représentant Autorisé dans la Communauté Européenne<br>Autorisierter Repräsentant in der EU | <b>IVD</b><br>In vitro diagnostic medical device<br>Dispositif médical de diagnostic in vitro<br>In-vitro-Diagnostikum  |  |

Dako North America, Inc.  
6392 Via Real  
Carpinteria, California 93013 USA

Tel 805 566 6655  
Fax 805 566 6688  
Technical Support 800 424 0021  
Customer Service 800 235 5763

**EC REP**

Dako Denmark A/S  
Produktionsvej 42  
DK-2600 Glostrup Denmark

Tel +45 4485 9500  
Fax +45 4485 9595

www.dako.com

PT0039/Rev C