**Monoclonal Mouse
Anti-Human
Calponin
Clone CALP**

Code M3556

Intended use

For In Vitro Diagnostic Use.

Summary and explanation

Calponin is a calmodulin, F-actin and tropomyosin binding protein which is thought to be involved in the regulation of smooth muscle contraction.^{2,3} Calponin expression is restricted to smooth muscle cells and has been shown to be a marker of the differentiated (contractile) phenotype of developing smooth muscle.²⁻⁵

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Anti-human calponin, CALP is a mouse monoclonal antibody supplied in liquid form as tissue culture supernatant (containing fetal bovine serum) dialyzed against 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide.

Clone: CALP¹ Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

Anti-calponin, CALP may be used at a dilution of 1:50 in the LSAB method determined on enzyme digested formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

Immunogen

Crude human uterus extract¹

Specificity

Multiple isoelectric variants of calponin have been identified, however only two molecular weight isoforms exist; a 34 kD form and a 29 kD form.⁵ Expression of the 29 kD form, *I*-calponin, is primarily restricted to muscle of the urogenital tract, whereas the higher molecular weight variant has been demonstrated in vascular and visceral smooth muscle.²⁻⁵ In Western blotting, monoclonal anti-calponin, CALP, has been demonstrated to react with only the 34 kD form of calponin in extracts of human aortic medial smooth muscle and is unreactive with fibroblast extracts of cultivated human foreskin.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-calponin, CALP can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes should be performed prior to staining.

For improved staining results, the deparaffinized tissue sections may be treated with a proteolytic enzyme followed by target retrieval. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended. Deparaffinized tissue sections should first be treated for 5 to 10 minutes with a mild enzyme solution. A recommended proteolytic enzyme is Proteinase K (code S3004) which must be further diluted 1:500 in 0.05 mol/L Tris-HCl, pH 7.6 to give a final concentration of 0.04 mg/mL.

Following proteolytic digestion, tissue sections can be treated with heat. When using the water bath method, preheat a Coplin jar containing 0.01 mol/L citrate buffer, pH 6.0 as well as a water bath to 95–99 °C. When the temperature has stabilized, place tissue sections into the Coplin jar containing the preheated buffer. Heat the tissue sections for 40 minutes. For improved staining results and a shorter incubation time, Target Retrieval Solution (code S1700) can be used in place of the 0.01 mol/L citrate buffer. Under these conditions the incubation time in the water bath may be reduced to 20 minutes. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with distilled water and place slides into buffer.

Cryostat Sections and Cell Smears

Anti-calponin, CALP can also be used to label cryostat sections or cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The staining pattern for anti-calponin is cytoplasmic.

Performance characteristics

Normal Cells

In immunohistochemical (IHC) studies on cryostat sections of human fetal tissue, monoclonal anti-calponin, CALP, was found to stain developing visceral smooth muscle of trachea, jejunum, esophagus and uterus in 10 and 20 week-old fetuses. Monoclonal anti-calponin did not react positively with 10 and 20 week-old fetal aortic smooth muscle cells.¹ Monoclonal antibody CALP was found to localize calponin in cryostat sections of adult visceral and vascular smooth muscle but not in epithelial cells, endothelial cells, or connective tissue fibroblasts. Adult aortic cells of the tunica media and a portion of subendothelial intimal cells were found to stain positively.¹ In cryostat sections and routinely fixed specimens of normal human breast, calponin expression has been demonstrated in smooth muscle cells of blood vessels and myoepithelial cells in the lobules, ducts and galactophorous sinuses.^{6,7} Periacinar and periductal myoepithelial cells of the salivary gland have also been shown to react positively with anti-calponin, whereas ductal epithelial cells were negative.⁸

Tumor Cells

Calponin expression has been demonstrated by IHC in myoepithelial cells in benign and malignant breast lesions.^{6,7} Myoepithelial cells in papillomas of the breast were found to stain positively with anti-calponin but no reactivity was observed in intracystic papillary carcinomas. In ductal carcinoma *in situ*, calponin immuno-reactivity has been demonstrated in myoepithelial cells at the periphery of involved ducts and lobules in complex sclerosing lesions of the breast.⁷ Anti-calponin was shown to label stromal myofibroblasts in the majority of invasive breast carcinomas tested but was unreactive with tumor cells.^{6,7} Reactivity has also been observed in neoplastic myoepithelium of routinely fixed pleomorphic adenomas of the salivary gland.⁸

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Références

Literatur

1. Frid MG, et al. Phenotypic changes of human smooth muscle cells during development: Late expression of heavy caldesmon and calponin. *Dev Biol* 1992; 153:185
2. Takahashi K, et al. Isolation and characterization of a 34,000-dalton calmodulin- and F-actin-binding protein from chicken gizzard smooth muscle. *Biochem Biophys Res Commun* 1986; 141(1):20
3. Takahashi K, et al. Occurrence of anti-gizzard p34K antibody cross-reactive components in bovine smooth muscles and non-smooth muscle tissues. *Life Sci* 1987; 41:291
4. Gimona M, et al. Smooth muscle specific expression of calponin. *FEBS Lett* 1990; 274(1,2):159
5. Draeger A, et al. Calponin: Developmental isoforms and a low molecular weight variant. *FEBS Lett* 1991; 291(1):24
6. Lazard D, et al. Expression of smooth muscle-specific proteins in myoepithelium and stromal myofibroblasts of normal and malignant human breast tissue. *Proc Natl Acad Sci USA* 1993; 90:999
7. Wang NP, et al. Monoclonal antibodies (MAbs) to novel myoepithelium-associated proteins can distinguish between benign and malignant lesions of the breast. *US Canad Acad Pathol Abstr* 1996; 26A:135
8. Savera AT, et al. Immunolocalization of smooth muscle-specific proteins (SMSPs) in pleomorphic adenomas. *US Canad Acad Pathol Abstr* 1996; 104A:603

REF Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use <i>Consulter les instructions d'utilisation</i> <i>Gebrauchsanweisung beachten</i>
 Manufacturer Fabricant Hersteller	LOT Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis
EC REP Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	IVD In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

EC REP

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0039/Rev C