**Monoclonal Mouse
Anti-Human
Estrogen Receptor α
Clone 1D5
Code M7047**

ENGLISH	
Intended use	For in vitro diagnostic use Monoclonal Mouse Anti-Human Estrogen Receptor α , Clone 1D5, is intended for use in immunohistochemistry. The antibody labels estrogen receptor α -positive cells and is useful in the assessment of estrogen receptor α status in human breast carcinomas (1-3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Steroid hormone receptors exhibit a high affinity and specificity for their ligands. The human estrogen receptor (ER) exists in two forms, ER α and ER β . The receptors have a high degree of homology in their DNA binding and ligand binding domains and can form heterodimers. However, the two ERs signal in different ways; for example, ER β does not potentiate the nuclear transcription factor AP-1 when complexed with oestradiol, whereas ER α does. This suggests that ER α and ER β may have different roles in gene regulation and their relative levels within tissues may influence cellular responses to oestrogens (3, 4). In breast cancer, the measurement of ER α status is generally recognized as being useful because the presence of ER α is related to an improved overall survival and a favourable response to endocrine therapy (2, 3, 5, 6).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> 1D5 (7). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial. The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
Immunogen	Soluble recombinant human estrogen receptor (7).
Specificity	The antibody specifically reacts with ER α and shows no reaction with ER β (8, 9). The epitope is located in the N-terminal domain (A/B region) of ER α . In Western immunoblotting using protein extracts from COS cells transfected with a plasmid vector expressing ER, the human breast cancer cell line MCF-7, and human endometrial tissue, the antibody labels a polypeptide of approximately 67 kDa corresponding to ER. No labelling of COS cells transfected with plasmid containing a deleted mutant with missing sequences encoding for the A/B region of ER was observed (7). As demonstrated by immunocytochemistry, the antibody reacts with the ER α -equivalent protein in rat (8).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is required. Optimal results are obtained with Dako Target Retrieval Solution, pH 9.0, code S2368, or Dako Target Retrieval Solution, High pH, code S3308. Less optimal results are obtained with Dako Target Retrieval Solution, code S1700, or 10 mmol/L citrate, pH 6.0. Pre-treatment of tissues with Dako Proteinase K, code S3020, was found inefficient. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling frozen sections (7).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human Estrogen Receptor α , code M7047, may be used at a dilution of 1:35 when applied on formalin-fixed, paraffin-embedded sections of human breast carcinoma and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval Solution, pH 9.0, code S2368, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, code S0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> Dako LSAB™+/HRP kit, code K0679, and Dako EnVision™+/HRP kits, codes K4004 and K4006, are recommended for frozen sections and cell preparations, the Dako APAAP kit, code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the Dako Autostainer
Staining interpretation	Cells labelled by the antibody display a nuclear staining pattern. Cytoplasmic labelling, if observed, should be regarded non-specific. In relation to tissues from human breast carcinomas, a variety of staining scales for ER-positivity has been reported in the literature. For example, staining results have been reported using either the percentage of positive cells or a staining intensity scale, or using the staining intensity scale in conjunction with the percentage of positive cells (H score) (2, 5, 6). A detailed description of different scoring systems is found in (6).

(104482-003)

M7047/EFG/KRM/2008.09.30 p. 1/4

Product specific limitations	1. Wash buffers containing high levels of detergent can decrease the staining intensity. 2. Occasional lymphoid tumours and non-lymphoid neoplasms such as melanomas are labelled.
Performance characteristics	<u>Normal tissues:</u> The antibody has been tested on a large range of normal tissues. Positive labelling was observed in the mammary gland, cervix uteri (squamous epithelium), uterus (endometrium) and prostate (cytoplasm of the fibromuscular stroma). <u>Abnormal tissues:</u> Numerous studies on formalin-fixed, paraffin-embedded breast cancer tissue sections have shown Anti-ER α , Clone 1D5, to be reliable and effective for the demonstration of ER α status (1, 2, 5, 6). The antibody also provides accurate prognostic information regarding response to endocrine therapies (2, 5, 6). Sensitivity data of 90% (2) and 89% (5), and specificity data of 51% (2) and 73% (5) have been reported when antibody labelling was compared with response to tamoxifen therapy. In (2), an H-score method with an arbitrary cut-off for positivity of 50 was used, while in (5) a cut-off point of 10% positive cells was used.

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