**Monoclonal Mouse
Anti-Human
CD235a, Glycophorin A**
Clone JC159
Code No./ Code/ Code-Nr. M 0819
Edition/ Ausgabe 27.02.03

ENGLISH	
Intended use	For in vitro diagnostic use. Monoclonal Mouse Anti-Human CD235a, Glycophorin A, Clone JC159, is intended for use in immunocytochemistry. The antibody labels erythroid cells at almost all stages of differentiation and is a useful tool for the identification of erythroleukaemia (1). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonyms for antigen	CD235a, GPA, Sialoglycoprotein alpha (1, 2).
Introduction	Glycophorin A is a single-pass membrane sialoglycoprotein and carrier of the blood group specificities M and N. The gene is located on chromosome 4q28-q31 and share sequence homology with genes encoding glycophorin B and glycophorin E (2). Glycophorin A is expressed by erythroid cells beginning on morphologically recognizable erythroid precursors, just after the CFU-E stage, to the mature erythrocyte. Once glycophorin A is maximally expressed, the amount in each erythroid cell remains constant and does not appear to change during further maturation (3). The majority of cases of erythroleukaemia express glycophorin A on neoplastic erythroblasts, whereas acute myeloid leukaemia and acute lymphoblastic leukaemia only very rarely express glycophorin A (1, 4).
Reagent provided	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN ₃ . <u>Clone:</u> JC159 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> See label on vial.
Immunogen	Membrane preparations from splenic hairy cell leukaemia cells (1).
Specificity	As demonstrated by haemagglutination and immunoblotting, the antibody reacts with a formalin-fixed, paraffin-embedded resistant epitope on the extracellular domain of the protein, probably located between amino acids 27 and 40. The antibody strongly labels normal erythroid cells at all stages of differentiation from the erythroblast to the mature red cell and does not react with glycophorin B (1).
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
Specimen preparation	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with proteinase K or heat-induced epitope retrieval is recommended. For heat-induced epitope retrieval optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed, frozen sections and cell smears (1).
Staining procedure	<u>Dilution:</u> Monoclonal Mouse Anti-Human CD235a, Glycophorin A, code No. M 0819, may be used at a dilution range of 1:50-1:400 when applied on formalin-fixed, paraffin-embedded sections of human tonsil and using 20 minutes heat-induced epitope retrieval in 10 mmol/L citrate buffer, pH 6.0, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. For frozen sections and cell preparations, the DakoCytomation APAAP kit, code No.
(102464-001)	M 0819/EFG/COS/27.02.03 p. 1/4

K 0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit.

Automation: The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.

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