**Monoclonal Mouse
Anti-Human
Laminin
Clone 4C7**

ENGLISH
Code M0638

Intended use

For In Vitro Diagnostic Use.

Monoclonal Mouse Anti-Human Laminin, clone 4C7 is intended for laboratory use to identify qualitatively by light microscopy laminin positive cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and explanation

The extracellular matrix (ECM) is a complex three-dimensional scaffold which surrounds and separates most tissue compartments, determines tissue architecture and dynamically interacts with its surrounding cells. The basement membrane is a part of the ECM and is composed of collagens, proteoglycans, elastin and several glycoproteins of which laminin is the most abundant. Laminins are a family of multifunctional basement membrane glycoproteins with many isoforms composed of variant alpha, beta, and gamma chains.² The laminins consist of at least 11 polypeptides (5 alpha chains, 3 beta chains, and 3 gamma chains) specific to basement membranes.³

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

M0638

Anti-Laminin, 4C7 is a monoclonal mouse antibody supplied in liquid form as ascites in 0.05 mol/L Tris-HCl, pH 7.6, 0.015 mol/L sodium azide and stabilizing protein.

Clone: 4C7 Isotype: IgG_{2a}, kappa
Mouse IgG concentration mg/L: See label on vial.

Anti-Laminin, 4C7 can be used at a dilution of 1:20 to 1:25 in the LSAB method determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

Immunogen

human placental laminin¹

Specificity

Anti-laminin, 4C7 reacts with a 380 kD laminin alpha5 chain (this antibody was previously thought to react with a 400 kD laminin alpha1 chain).^{1,5} Anti-laminin, 4C7 was shown to partially block the neurite-promoting activity of laminin.^{4,6}

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation

Paraffin Sections

Anti-Laminin, 4C7 can be used on formalin-fixed, paraffin-embedded tissue sections. Pretreatment of tissue with proteolytic enzymes should be performed before staining.

Cryostat Sections

Anti-Laminin, 4C7 can also be used to label laminin in cryostat sections.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-laminin is cytoplasmic and/or membranous.

Performance characteristics

Normal Tissues

Anti-laminin, 4C7 stained proximal and distal tubuli, glomerular basement membrane, and blood vessels of adult kidney; basement membrane of normal muscle and blood vessels of muscle tissue; and placenta.¹ Anti-laminin, 4C7 stained muscle fibers in addition to blood vessels in human fetal muscle. 4C7 antigen is widely distributed in developing and adult human, rat, and hamster tissues.¹

Abnormal Tissues

Anti-laminin, 4C7 demonstrated increased staining in muscular dystrophies (basement membranes of muscle fibers, blood vessels, and somewhat more intensely, groups of small diameter regenerating muscle fibers).¹

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Bibliographie

Literaturangaben

1. Tiger CF, Champlaud MF, Pedrosa-Domellof F, Thornell LE, Ekblom P, Gullberg D. Presence of laminin alpha5 chain and lack of laminin alpha1 chain during human muscle development and in muscular dystrophies. *J Biol Chem* 1997; 272(45):28590
2. Jaspars LH, De Melker AA, Bonnet P, Sonnenberg A, Meijer CJ. Distribution of laminin variants and their integrin receptors in human secondary lymphoid tissue. Colocalization suggests that the alpha 6 beta 4-integrin is a receptor for laminin-5 in lymphoid follicles. *Cell Adhes Commun* 1996; 4(4–5):269
3. Makino M, Okazaki I, Kasai S, Nishi N, Bougaeva M, Weeks BS, Otaka A, Nielsen PK, Yamada Y, Nomizu M. Identification of cell binding sites in the laminin alpha5-chain G domain. *Exp Cell Res* 2002; 277(1):95
4. Engvall E, Earwicker D, Haaparanta T, Ruoslahti E, Sanes JR. Distribution and isolation of four laminin variants; tissue restricted distribution of heterodimers assembled from five different subunits. *Cell Regulation* 1990; 1:731
5. Church HJ, Aplin JD. BeWo choriocarcinoma cells produce laminin 10. *Biochem J*. 1998; 332(pt2):491
6. Kleinman HK, Ogler RC, Cannon FB, Little CD, Sweeney TM, Luckenbill-Edds L. Laminin receptors for neurite formation. *Proc Natl Acad Sci USA* 1988; 85:1282

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Contains sufficient for <N> tests Contenu suffisant pour <n> tests Inhalt ausreichend für „n“ Ansätze
 Use by Utiliser jusque Verwendbar bis	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU

Dako North America, Inc.
6392 Via Real
Carpinteria, California 93013 USA

Tel 805 566 6655
Fax 805 566 6688
Technical Support 800 424 0021
Customer Service 800 235 5763

Dako Denmark A/S
Produktionsvej 42
DK-2600 Glostrup Denmark

Tel +45 4485 9500
Fax +45 4485 9595

www.dako.com

PT0020/ Rev C