**Monoclonal Mouse  
Anti-Human  
Ki-67 Antigen**  
Clone MIB-1  
Code M7240

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use.  Monoclonal Mouse Anti-Human Ki-67 Antigen, Clone MIB-1, is intended for use in immunohistochemistry. With more than 400 literature citations, the MIB-1 antibody has now been established as a reference monoclonal mouse antibody for the demonstration of the Ki-67 antigen in formalin-fixed, paraffin-embedded specimens. In diagnostic histopathology and cell biology, the antibody has proven valuable for the demonstration of the Ki-67 antigen in normal and neoplastic cells, for example in soft-tissue sarcoma (1), prostatic adenocarcinoma (2), and breast carcinoma (3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Introduction</b>	The Ki-67 antigen is a nuclear protein, which is defined by its reactivity with monoclonal antibody from the Ki-67 clone (4). Two isoforms of 345 and 395 kDa have been identified (5). The Ki-67 antigen is preferentially expressed during all active phases of the cell cycle (G <sub>1</sub> , S, G <sub>2</sub> and M-phases), but it is absent in resting cells (G <sub>0</sub> -phase) (4). During interphase, the antigen can be exclusively detected within the nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. The antigen is rapidly degraded as the cell enters the non-proliferative state (6), and there appears to be no expression of Ki-67 during DNA repair processes (7).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, 1% Bovine Serum Albumin, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> MIB-1 (8). <u>Isotype:</u> IgG1, kappa. <u>Mouse Ig concentration:</u> see label on vial.  The protein concentration between lots may vary without influencing the optimal dilution. The titer of each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.
<b>Immunogen</b>	Human recombinant peptide corresponding to a 1002 bp Ki-67 cDNA fragment (8).
<b>Specificity</b>	In Western blotting of lysates of the multiple myeloma cell line, IM-9, the MIB-1 antibody labels bands of 345 and 395 kDa, identical to the bands labelled by the original Ki-67 antibody. Furthermore, Western blotting and competitive binding experiments clearly demonstrate that MIB-1, like the original Ki-67 antibody, reacts with an epitope encoded by a 66 bp repetitive element in the Ki-67 gene. In immunohistochemistry, the MIB-1 and the Ki-67 antibodies provide identical staining patterns on serial tonsillar frozen sections (8). The MIB-1 antibody recognizes native Ki-67 antigen and recombinant fragments of the Ki-67 molecule (8).  As demonstrated by immunohistochemistry, the antibody cross-reacts with the Ki-67-equivalent protein in various mammals, including cow, dog, horse, sheep (6), and swine.
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labeling paraffin-embedded tissue sections fixed in formalin, Bouin's fixative, B5-fixative, Dubosq-Brazil or Zenker's solution (9). For tissues fixed in formalin, optimal results are obtained with Dako Target Retrieval Solution, Code S1700, Dako Target Retrieval Solution, High pH, Code S3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris/HCl, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found inefficient. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. <u>Frozen sections and cell preparations:</u> The antibody can be used for labelling acetone-fixed frozen sections and fixed cell smears.
<b>Staining procedures</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human Ki-67 Antigen, Code M7240, may be used at a dilution range of 1:75-1:150 when applied on formalin-fixed, paraffin-embedded sections of human tonsil or human intestinal mucosa and using 20 minutes heat-induced epitope retrieval in Dako Target Retrieval solution, Code S1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is Dako Mouse IgG1, Code X0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in Dako Antibody Diluent, Code S0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> Dako LSAB™+/HRP kit, Code K0679, and Dako EnVision™+/HRP kits, Codes K4004 and K4006, are recommended. For frozen sections and cell preparations, the Dako APAAP kit, Code K0670, is a good alternative if endogenous peroxidase staining is a concern. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunohistochemical staining using automated platforms, such as the Dako Autostainer.

(104980-003)

M7240/EFG/KRM/2008.11.07 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

<b>Product-specific limitations</b>	Occasional labelling of tissue components in vessel walls and pancreatic stroma has been observed in immunohistochemistry.
<b>Performance characteristics</b>	Cells labelled by the antibody display a nuclear staining pattern except in mitotic cells, where the chromosomes and the cytoplasm are labelled. <u>Normal tissues:</u> In small intestine and colon, mucosal neck cells and deep crypt cells are positive with the antibody. The same is the case with germinal centre cells of Payer's patches. Surface epithelium cells and other mucosal cells, e.g. Paneth's cells and cells of submucosal connective tissue are completely negative. The same is true for muscle cells, with the exception of a few positive cells in smooth muscles. Kidney, liver, pancreas and brain are negative (9). <u>Abnormal tissues:</u> In a study of 123 soft tissue sarcomas, Ki-67 was demonstrated with the highest expression in malignant fibrous histiocytomas and the lowest in liposarcomas using the MIB-1 antibody. The antibody was also applied successfully for the demonstration of the Ki-67 antigen in 221 prostate carcinomas (2) and 919 breast carcinomas (3).

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