**Monoclonal Mouse
Anti-MyoD1
Clone 5.8A**

ENGLISH
Code M3512**Intended use**

For In Vitro Diagnostic Use.

Monoclonal mouse anti-MyoD1 is intended for laboratory use to identify qualitatively by light microscopy MyoD1 positive cells in normal and neoplastic tissues using immunohistochemical (IHC) test methods. The clinical interpretation of any positive staining or its absence should be complemented by morphological and histological studies with proper controls. Evaluations should be made within the context of the patient's clinical history and other diagnostic tests by a qualified individual.

Summary and explanation

The MyoD1 gene belongs to a family of myogenic determination genes, including myogenin, myf-5 and MRF4, all of which encode transcription factors.^{2,3} Transfection of MyoD1 cDNA into non-muscle cells has been shown to activate expression of muscle specific genes and in some cases induce myogenesis.^{4,5} The MyoD1 protein is a 45 kD nuclear phosphoprotein which induces myogenesis through transcriptional activation of muscle-specific genes.^{3,6,7} Nuclear expression of MyoD1 is restricted to skeletal muscle tissue and has been demonstrated to be a sensitive marker of myogenic differentiation.^{1,8,9}

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Monoclonal Mouse antibody provided in liquid form as tissue culture supernatant in 0.05 mol/L Tris-HCl, pH 7.2 and 0.015 mol/L sodium azide. This product contains stabilizing protein.

Clone: 5.8A Isotype: IgG₁, kappa
Mouse IgG concentration mg/L: See label on vial.

Monoclonal anti-MyoD1 may be used at a dilution of 1:50 to 1:75 in the LSAB™+ HRP method, determined on formalin-fixed, paraffin-embedded tissue. These are guidelines only; optimal dilutions should be determined by the individual laboratory.

ImmunogenRecombinant wild-type murine MyoD1 protein¹**Specificity**

Monoclonal anti-MyoD1 has been shown by Western blot and immunoprecipitation assays to bind an epitope between amino acids 180–189 on the carboxy terminus of the MyoD1 protein.¹

Materials required, but not supplied

Refer to Dako's *General Instructions for Immunohistochemical Staining* and/or the detection system instructions.

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, NaN₃ may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused reagents should be disposed of according to local, State, and Federal regulations.

Storage

Store at 2–8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation**Paraffin Sections**

Anti-MyoD1 can be used on formalin-fixed, paraffin-embedded tissue sections.

The deparaffinized tissue sections must be treated with heat prior to the immunohistochemical (IHC) staining procedure. For greater adherence of tissue sections to glass slides, the use of Silanized Slides (code S3003) is recommended.

When using the water bath method, preheat a Coplin jar containing 0.01 mol/L citrate buffer, pH 6.0 as well as a water bath to 95–99 °C. When the temperature has stabilized, place tissue sections into the Coplin jar containing the preheated buffer. Heat the tissue sections for 40 minutes. For improved staining results, Target Retrieval Solution (code S1700) can be used in place of the 0.01 mol/L citrate buffer. After thermal treatment, allow the jar with buffer and slides to cool for 20 minutes at room temperature. Rinse well with distilled water and place slides into buffer.

For this antibody to perform optimally on paraffin-embedded tissues, a high sensitivity detection system is recommended such as the LSAB™+ (code K0679).

Cryostat Sections and Cell Smears

Anti-MyoD1 can also be used to label cryostat sections or cell smears.

Staining procedure

Follow the recommended procedure for the detection system selected.

Staining interpretation

The cellular staining pattern for anti-MyoD1 is nuclear.

Performance characteristics

Normal Cells

In immunohistochemical assays, monoclonal anti-MyoD1 strongly stains the nuclei of myoblasts in developing skeletal muscle tissue whereas the majority of adult skeletal muscle has been found to be negative.^{1,8} Monoclonal anti-MyoD1 has been tested in immunohistochemistry on a variety of additional normal tissues and was not shown to label the nuclei of any of the following tissues: adrenal gland, brain, breast, fetal skeletal muscle (16 weeks), heart, liver, lung, lymph node, ovary, pancreas, skin, small intestine, spleen, thymus, and thyroid.¹ Weak cytoplasmic staining has been observed in several non-muscle tissues, including glandular epithelium.⁸

Tumor Cells

MyoD1 immunostaining has been demonstrated in the majority of rhabdomyosarcomas of various histological subtypes.^{1,8} Results of a recent study by Wang et al suggest that MyoD1 expression in rhabdomyosarcomas is inversely related to the degree of cellular differentiation of the tumor cells.⁸ Nuclear expression of MyoD1 has been demonstrated in other tumors, including ectomesenchymoma and a subset of Wilm's tumors.¹ Cytoplasmic immunoreactivity with monoclonal anti-MyoD1 has been reported in rhabdomyosarcomas, neuroblastomas, Ewing's sarcomas and alveolar soft part sarcomas. However, only nuclear staining should be considered as evidence of skeletal myogenic differentiation.^{8,10}

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Références

Literatur

1. Dias P, Parham DM, Shapiro DN, Tapscott SJ, Houghton PJ. Monoclonal antibodies to the myogenic regulatory protein MyoD1: Epitope mapping and diagnostic utility. *Canc Res* 1992;52(23):6431-9
2. Weintraub H, Davis R, Tapscott S, Thayer M, Krause M, Benezra R, Blackwell TK, Turner D, Rupp R, Hollenberg S, Zhuang Y, Lassar A. The MyoD gene family: Nodal point during specification of the muscle cell lineage. *Science* 1991;251(4995):761-6
3. Dias P, Dilling M, Houghton P. The molecular basis of skeletal muscle differentiation. *Seminars in Diagnostic Pathology* 1994;11(1):3-14
4. Davis RL, Weintraub H, Lassar AB. Expression of a single transfected cDNA converts fibroblasts to myoblasts. *Cell* 1987;51(6):967-1000
5. Weintraub H, Tapscott SJ, Davis RL, Thayer MJ, Adam MA, Lassar AB, Miller AD. Activation of muscle-specific genes in pigment, nerve, fat, liver, and fibroblast cell lines by forced expression of MyoD. *Proc Natl Acad Sci USA* 1989;86(14):5434-8
6. Davis RL, Cheng PF, Lassar AB, Weintraub H. The MyoD DNA binding domain contains a recognition code for muscle-specific gene activation. *Cell* 1990;60(5):733-46
7. Tapscott SJ, Davis RL, Thayer MJ, Cheng PR, Weintraub H, Lassar AB. MyoD1: A nuclear phosphoprotein requiring a Myc homology region to convert fibroblasts to myoblasts. *Science* 1988;242(4877):405-11
8. Wang NP, Marx J, McNutt MA, Rutledge JC, Gown AM. Expression of myogenic regulatory proteins (myogenin and MyoD1) in small blue round cell tumors of childhood. *Amer J Path* 1995;147(6):1799-810
9. Tallini G, Parham DM, Dias P, Cordon-Cardo C, Houghton PJ, Rosai J. Myogenic regulatory protein expression in adult soft tissue sarcomas. *Amer J Path* 1994;144(4):693-701
10. Wang NP, McNutt MA, Bacchi CE, Gown AM. Expression of myogenic regulatory proteins (myogenin and MyoD1) and muscle-specific cytoskeleton proteins in alveolar soft part sarcoma. *Mod Path* 1995;8(1):13A

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA Tel 805 566 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763	 Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel +45 4485 9500 Fax +45 4485 9595 www.dako.com
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Use by Utiliser jusque Verwendbar bis		
 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum			

PT0039/Rev C

Edition 04/07