**Monoclonal Mouse Anti-Human Neuron-Specific Enolase**  
Clone BBS/NC/VI-H14  
**Code No./ Code/ Code-Nr. M 0873**  
Edition/ Ausgabe 20.01.03

ENGLISH	
<b>Intended use</b>	For in vitro diagnostic use. Monoclonal Mouse Anti-Human Neuron-Specific Enolase (NSE), Clone BBS/NC/VI-H14, is intended for use in immunocytochemistry. The antibody labels both normal and neoplastic cells of neuronal and neuroendocrine origin (1), and although NSE is not an exclusive neuronal marker, it may be used for the identification of peripheral nerves, neural and neuroendocrine tumours, such as neuroblastomas, retinoblastomas, desmoplastic malignant melanoma and small-cell lung cancer (2, 3). Differential identification is aided by the results from a panel of antibodies. Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
<b>Synonyms for antigen</b>	$\gamma$ -enolase and 14-3-2 protein (4)
<b>Introduction</b>	Enolases (2-phospho-D-glycerate hydrolase) are glycolytic homo- or heterodimeric isoenzymes, which catalyse the interconversion of 2-phosphoglycerate and phosphoenolpyruvate. Three subunits, $\alpha$ (46 kDa), $\beta$ (44 kDa) and $\gamma$ (46 kDa), constitute the building blocks of the enolase molecules, and five forms ( $\alpha\alpha$ , $\beta\beta$ , $\gamma\gamma$ , $\alpha\beta$ and $\alpha\gamma$ ) of the enzyme have been identified (5). The enolases are divided in three main groups according to their localization: non-neuronal enolase (NNE or $\alpha$ ), muscle-specific enolase (MSE or $\beta$ ), and neuron-specific enolase (NSE or $\gamma$ ) (3).
<b>Reagent provided</b>	Monoclonal mouse antibody provided in liquid form as cell culture supernatant dialysed against 0.05 mol/L Tris/HCl, pH 7.2, and containing 15 mmol/L NaN <sub>3</sub> . <u>Clone:</u> BBS/NC/VI-H14 (1). <u>Isotype:</u> IgG1, kappa. <u>Mouse IgG concentration:</u> see label on vial.
<b>Immunogen</b>	$\gamma\gamma$ -enolase purified from human brain by a dissociation-reassociation method in a chromatofocusing column (1).
<b>Specificity</b>	In Western blotting of reduced, SDS-denatured purified human $\gamma\gamma$ -enolase, $\alpha\alpha$ -enolase, and extract of human muscle containing $\beta\beta$ -enolase, the antibody labels solely the $\gamma$ -enolase subunit (1). As demonstrated by immunocytochemistry, the antibody cross-reacts with the NSE-equivalent protein in guinea pig (1).
<b>Precautions</b>	1. For professional users. 2. This product contains sodium azide (NaN <sub>3</sub> ), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. This product contains material of animal origin, therefore the product should be handled as potentially infectious.
<b>Storage</b>	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.
<b>Specimen preparation</b>	<u>Paraffin sections:</u> The antibody can be used for labelling paraffin-embedded tissue sections fixed in formalin. Pre-treatment of tissues with heat-induced epitope retrieval is recommended. Optimal results are obtained with DakoCytomation Target Retrieval Solution, code No. S 1700, DakoCytomation Target Retrieval Solution, High pH, code No. S 3308, 10 mmol/L citrate buffer, pH 6.0, or 10 mmol/L Tris buffer, 1 mmol/L EDTA, pH 9.0. Pre-treatment of tissues with proteinase K was found destructive of the epitope. The tissue sections should not dry out during the treatment or during the following immunocytochemical staining procedure.
<b>Staining procedure</b>	<u>Dilution:</u> Monoclonal Mouse Anti-Human Neuron-Specific Enolase (NSE), code No. M 0873, may be used at a dilution range of 1:100-1:200 when applied on formalin-fixed, paraffin-embedded sections of human colon and using 20 minutes heat-induced epitope retrieval in DakoCytomation Target Retrieval Solution, code No. S 1700, and 30 minutes incubation at room temperature with the primary antibody. Optimal conditions may vary depending on specimen and preparation method, and should be determined by each individual laboratory. The recommended negative control is DakoCytomation Mouse IgG1, code No. X 0931, diluted to the same mouse IgG concentration as the primary antibody. Unless the stability of the diluted antibody and negative control has been established in the actual staining procedure, it is recommended to dilute these reagents immediately before use, or dilute in DakoCytomation Antibody Diluent, code No. S 0809. Positive and negative controls should be run simultaneously with patient specimen. <u>Visualization:</u> DAKO LSAB™+/HRP kit, code No. K 0679, and DAKO EnVision™+/HRP kits, code Nos. K 4004 and K 4006, are recommended. Follow the procedure enclosed with the selected visualization kit. <u>Automation:</u> The antibody is well-suited for immunocytochemical staining using automated platforms, such as the DakoCytomation Autostainer.
<b>Performance characteristics</b>	Cells labelled by the antibody display a cytoplasmic staining pattern. In neurons, labelling occurs in both cytoplasm and processes.
(102353-001)	M 0873/EFG/HEW/20.01.03 p. 1/4

Normal tissues: The antibody labels neurons in human cerebral cortex and brainstem nuclei (1). Additionally, the antibody labels ganglion cells in the human gastrointestinal tract and myelinated and unmyelinated nerve fibers (1). In normal adult testis, the antibody displays weak immunostaining of the spermatogonia (3).

Abnormal tissues: The antibody labels tumours derived from or containing neurons, ganglion cells and peripheral nerves. However, also non-neuronal derived tumours, as meningiomas, medulloblastomas, astrocytomas, glioblastomas, oligoastrocytomas, oligodendrogliomas, pituitary adenomas, schwannomas, ependymomas, meningiosarcomas, gliosarcomas and metastases of different origin – i.e. melanomas, small-cell carcinomas of lung, and undifferentiated carcinomas may be labelled specifically (1, 4). In addition, the antibody labels germ cell tumours, i.e. testicular carcinomas-in-situ, seminomas and embryonal carcinomas (3).

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